

Golden Compass™ Sport (GP605).



Available accessories:

- Crutch/cane holder tube
- Oxygen tank holder
- Quad cane holder
- Walker holder
- Swing-away desk arm
- Swing-away footrests
- Stump support
- Elevating leg rests
- Power chair cover
- Rear basket
- Pack n' go
- Van/pan seat
- EZ lock system
- Universal cup holder
- Armrest bag

Golden Compass™ Sport specifications:

Weight capacity	300 lbs.
Maximum speed	4 mph
Length	39"
Width	24"
Turning radius	19.5"
Range*	15 miles

Seat sizes:

Maximum seat dimensions	18" x 18"/20"
-------------------------	---------------

Technical:

Drive wheels	Center
Front casters	6"
Rear casters	6"
Ground clearance	3.25"
Drivetrain	13 months
Battery	(2) U1 batteries
Battery charger	Off-board
Battery weight	22 lbs. each

Colors:	Red	Blue
----------------	-----	------

The Golden Compass Sport features the new Dynamic Linx controller! The new Linx offers a better ride experience using patented Dynamic Load Compensation technology. This new technology improves the terrain-climbing ability of the Compass Sport's patented springless articulating chassis. Learning to drive the Compass Sport is easier and a more natural experience than before!

*Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

To pre-qualify for a power chair or scooter at little or no cost to you.
 Call toll-free: 866.387.2668, lincare.com

