

Pride Jazzy® 1450.



Features:

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Oversized rear casters for optimal performance and comfort
- Universal, fixed, in-line joystick mount
- Large 17" x 13" foot platform
- Superb traction, flat-free, drive tires
- Extremely efficient motor package

Pride Jazzy® 1450 specifications:

Weight capacity	600 lbs.
Maximum speed ¹	Up to 4.5 mph
Length ²	38.25"
Width ²	29.25"
Turning radius	26.5"
Range ^{2,4}	Up to 15.7 miles

Seat sizes:

Standard seat size	24" x 24"
--------------------	-----------

Technical:

Drive wheels	14" solid
Rear casters	9"
Ground clearance ²	1.5"
Drivetrain	Two-motor, front-wheel drive
Battery ^{3,5}	(2) 12V, group 24
Battery charger	Off-board, 8 amp
Battery weight	52 lbs. each

Colors:

Candy Apple Red	Viper Blue
-----------------	------------

The Jazzy® 1450 offers heavy-duty construction and performance components. This bariatric power chair features a 600 lb. weight capacity, front-wheel drive for superb indoor maneuverability in the tightest of spaces and large front drive wheels for extra absorption over obstacles. The Jazzy 1450 combines superb climbing capabilities with front anti-tips for transitioning over grade changes and handling various terrains easily.

To pre-qualify for a power chair or scooter at little or no cost to you.
Call toll-free: 866.387.2668, lincare.com



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.