

# Pride Jazzy® Elite HD.



## Features:

- Depth-adjustable, limited-recline, fold-flat, highback seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 70" lap belt

## Pride Jazzy® Elite HD specifications:

Weight capacity	450 lbs.
Maximum speed <sup>1</sup>	4 mph
Length <sup>2</sup>	35.5"
Width <sup>2</sup>	24.25"
Turning radius	24"
Range <sup>2,4</sup>	Up to 9.8 miles

## Seat sizes:

High-back seat sizes	18" x 18-20", 20" x 18-20" and 22" x 20-22"
----------------------	---

## Technical:

Drive wheels	14" solid
Rear casters	6" solid casters
Front wheels	3" solid anti-tips
Ground clearance <sup>2</sup>	3"
Drivetrain	Two-motor, front-wheel drive
Battery <sup>3,5</sup>	(2) 12V, U-1
Battery charger	Off-board, 5 amp
Battery weight	37 lbs. each

## Colors:

Jazzy Red	Jazzy Blue
-----------	------------

The Jazzy® Elite HD features two-motor, front-wheel drive technology and heavy-duty construction. The Jazzy Elite HD features 14" drive wheels for enhanced outdoor performance combined with a frontwheel drive design to provide tight turns for optimal indoor maneuverability in small spaces. Large front wheels add extra absorption for superb climbing capabilities and combine with front anti-tip wheels for handling various terrains with ease.

To pre-qualify for a power chair or scooter at little or no cost to you.

Call toll-free: 866.387.2668, [lincare.com](http://lincare.com)



1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).  
 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).  
 3 AGM or gel cell type required.  
 4 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.  
 5 Battery weight may vary based on manufacturer.