

Pride® Jazzy® EVO 613 Series.



The most advanced Jazzy® ever. You will love the all-new Medicare-reimbursable Jazzy® EVO 613 series!

To pre-qualify for a power chair or scooter at little or no cost to you. Call toll-free: 866.387.2668, lincare.com

Features:

- Maximum speed of up to 5.1 mph
- 13" mid-wheel drive tires
- Narrow base width of 22"
- Group 2, standard (K0822/K0823)
- Choice of six contemporary removable, matte-finish color through shrouds
- Depth adjustable, limited-recline, memory foam, high back with single set post allows the user to swivel
- Pride's patented Active-Trac® (ATX) suspension
- Mid-Wheel 6® technology

Pride Jazzy® EVO 613, 613 Li (K0822, K0823) specifications:

Weight capacity	300 lbs.
Maximum speed ²	Up to 5.1 mph
Dimensions ³	Width: 22", Length: 41" w/ footrest
Total weight	EVO 613: 189 lbs
Turning radius	22"
Range per charge ³	EVO 613: Up to 14 miles/300 lbs. EVO 613: Up to 17.1 miles/200 lbs.

Warranty:

Frame	5-year limited
Drivetrain	2-year limited
Electronics	2-year limited
SLA battery	13-month limited

Matte-colors:



¹ Pride FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. These specifications can be subject to a variance of (+/- 10%).

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+/- 3%).