

Pride Jazzy® Select.



Features:

- Mid-wheel drive technology and a compact size is perfect for small spaces by offering a tight turning radius of 22.5"
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Articulating rear caster beam works in conjunction with the ATX Suspension for a smooth ride over uneven surfaces
- Easy-access free-wheel lever
- Color-through shroud keeps the shroud looking new by reducing the visual scrapes and chips

Pride Jazzy® Select specifications:

Weight capacity	300 lbs.
Maximum speed ¹	Up to 4.3 mph
Length ²	44"
Width ²	23.74"
Turning radius	22.5"
Range ^{2,4}	Up to 19.4 miles
Seat sizes:	
High-back seat sizes	16"x16"-18", 18"x18"-20", 20"x18"-20"
Technical:	
Drive wheels	10" solid
Rear wheels	6" solid
Front wheels	5" solid
Ground clearance ²	0.75" at motor bracket
Drivetrain	Two-motor, mid-wheel drive®
Battery ^{3,5}	(2) 12V, U-1
Battery charger	Off-board, 3.5 amp
Battery weight	25.2 lbs. each

Colors:

Jazzy Red

Experience the simplicity and maneuverability of the Jazzy Select®. Equipped with Active-Trac® Suspension and dual motors, this power chair ensures a smooth, comfortable ride over uneven terrain. Midwheel drive offers a tight turning radius, making the Jazzy Select a perfect choice both indoors and out.

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

To pre-qualify for a power chair or scooter at little or no cost to you.

Call toll-free: 866.387.2668, lincare.com